**BVDA INTERNATIONAL**

PO BOX 2323 2002 CH HAARLEM HOLLAND
TEL +31 (23) 5424708 FAX +31 (23) 5322358
EMAIL INFO@BVDA.NL WWW.BVDA.COM

Nastanak

I,2-Indandioni su prvo bili napravljeni kao agens za pripremu zamene nihidrina (npr. 5-metiltioninhidrina) od strane prof. Madeleine Joullié i njenih studenata postdiplomaca, Dr Diane Hauze i Dr Olga Petrovskaia, sa Univerziteta u Pensilvaniji (Upenn). Tokom provere sastojaka pri Tajnoj službi SAD koju su radili Dr Tony Cantu and Robert Ramotowski, otkriveno je da ove supstance imaju i svojstvo da latente otiske čine vidljivim.

Indandion u poređenju sa DFO

Nekoliko zamenjenih indandiona je napravljeno na Upenn i provereno kod T.SI. SAD [1]. Dr. Joseph Almog i istraživači iz Izraelske narodne policije su sintetisali i proverili znatne količine indandiona [2]. Australijski istraživači su uporedili I,2-indandion i 5,6-dimetoksi-I,2-indandion sa ninhidrinom i DFO [3]. Wiesner et al. su proverili performanse I,2-indandiona (2g/l u rastvoru HFE-7100 koji sadrži 7% etil acetat) kod korišćenih kontrolnih materijala [5].

Ovi tekstovi o istraživanju su imali za zaključak da indandioni mogu vizualizovati latentne otiske jednako dobro ili bolje od DFO. Kao i DFO, indandioni ne proizvode jako obojeni prvobitno razvijeni otisak (obično svetlo roze boje), ali daju jak fluoroscentan detalj kada se osvetle zelenim svetlom.

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
Comparison of the fluorescence of DFO (left) and 1,2-IND (green light/dark orange filter).

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Reference

- [1] Ramotowski, R.; Cantu, A.A.; Joullié, M.M.; Petrovskaia, O.; "1,2-Indanediones: A Preliminary Evaluation of a New Class of Amino Acid Visualizing Compounds," *Fingerprint World* **1997**, Vol. 23, No. 90, pp. 131-140.
- [2] Almog, J.; Springer, E.; Wiesner, S.; Frank, A. et al., "Latent Fingerprint Visualization by 1,2-Indanedione and Related Compounds: Preliminary Results," *Journal of Forensic Sciences* **1999**, Vol. 44, No. 1, pp. 114-118.
- [3] Roux, C.; Jones, N.; Lennard, C.; Stoilovic, M., "Evaluation of 1,2-Indanedione and 5,6-Dimethoxy-1,2-Indanedione for the Detection of Latent Fingerprints on Porous Surfaces", *Journal of Forensic Sciences* **2000**, Vol. 45, No. 4, pp. 761-769.
- [4] Unpublished observations by researchers from the University of Pennsylvania, USSS, Israel National Police, and PSDB.
- [5] Wiesner, S.; Springer, E.; Sasson, Y.; Almog, J., "Chemical Development of Latent Fingerprints: 1,2-Indanedione Has Come of Age", *Journal of Forensic Sciences* **2001**, Vol. 46, No. 5, pp. 1082-1084.
- [6] Merrick, S., Gardner, S.J.; Sears, V.G.; Hewlett, D.F., "An Operational Trial of Ozone-Friendly DFO and 1,2-Indanedione Formulations for Latent Fingerprint Detection", *Journal of Forensic Identification* **2002**, Vol. 52, No. 5, pp. 595-605.
- [7] Kasper, S.P.; Minnillo, D.J., Rockhold, A.M., "Validating IND (1,2 - indanedione)", *Forensic Science Communications* **2002**, Vol. 4, No. 4.
- [8] Stimac, J.T., "Thermal Paper: Latent Friction Ridge Development via 1,2-Indanedione", *Journal of Forensic Identification* **2003**, Vol. 53, No. 3, pp. 265-271.
- [9a] Wallace, Christie, "Characterisation and Optimisation of Reaction of 1,2-Indanedione with Amino Acids and its Application to Fingerprint Casework", presentation at the 17th International Symposium on the Forensic Sciences, March 2004 (abstract: <http://www.anzfss2004.org.nz/cgi-bin/ViewAbstracts.cgi?paperid=57>).
- [9b] Wallace, Christie; Claude, Roux; poster presented at the IAI-conference in St. Louis, August 2004 (to be published).